



# Smooth and Relax

The **hydraulic mechanism** ensures smooth and quiet movement of the joints and provides a comfortable feeling for patients.

## Features



Unique low patient entry seat position along with a **200 kg patient lifting capability** is one of the most versatile chairs in the market place.

## Color variations

The colors come in four themes: warm, cool, earth and vitamin. Any choices from these leather colors should match the interior appointments of your clinic. The upholstery can be chosen from seam, seamless and luxury type.



### [ Seam & Seamless ]

<b>Warm Color</b>						
	Coral Pink	Light Pink	Beige	Camel	Bordeaux	Black
<b>Cool Color</b>						
	Mint Blue	Indigo Blue	Light Blue	Purple	Light Gray	Gray Teal
<b>Earth Color</b>						
	Rose	Mauve	Forest	Aldo Blue	Ash Gray	Chocolate
<b>Vitamin Color</b>						
	Scarlet Red	Mandarin Orange	Lemon Yellow	Yellow Green	Sky Blue	Blue

### [ Luxury ]

Black	Dark Blue	Maroon	England Green	Orange	Camel
Mustard Yellow	Ivory				

## CLESTA II UNIT

### Compact unit supports versatile treatment styles

- Programmed micromotor setting ( E Type )
- Programmed scaler setting ( E Type )
- Doctor preset setting for 4 groups of combination
- Large operator's table
- Light: EURUS LIGHT  
900 LED LIGHT  
300 LED LIGHT

Unit : CLESTA II Holder Type  
Chair : CLESTA II ( EURUS Type )  
Light : EURUS LIGHT




### Double articulating headrest

The precise positioning of the headrest can be set in **7.5 degree increments** accommodating a wide range of patient types. The height adjustment for taller patients can be easily positioned with **very little effort** from the operator.



### Rotating armrest

The armrests are designed with a length and shape that makes it easy for patients to hold and the **side-rotating mechanism** helps elderly patients mount and dismount the chair with ease.



### Chair stick switch

The dental chair can be entirely manipulated in a **hands-free**, hygienic and intuitive manner using only the operator's foot ( patent pending ).